

THERMAL CONDUCTIVITY APPARATUS®

OBJECTIVE

Thermal Conductivity of Liquids & Nanofluids & Engine Oils etc.

THEORY

The Apparatus follows widely accepted theory of heat conduction in liquids based on Debye's concept in which the hydroacoustic vibrations (phonons) of a continuous medium (base fluid) are responsible for the heat transfer in liquids. Based on this heat transfer mechanism, Bridgman has obtained a formula, characterized by the direct proportionality between thermal conductivity and sound velocity in pure liquids.

$$k = 3.0 \left(\frac{N}{V} \right)^{\frac{2}{3}} K_B v$$

where v_s = ultrasound velocity,

N (Avogadro's number) = 6.02×10^{23} and V (molar volume) = m/ρ

K_B = (Boltzmann's constant) = 1.3807×10^{-23} J/K

For Polyatomic Liquids it was modified as:

$$k = 2.8 \left(\frac{N}{V} \right)^{\frac{2}{3}} k_B v$$

For Nanofluids it was further modified by J.Hemalatha as under:

$$k_{bm} = 2.8 \left(\frac{\rho_{nf} N_A}{M_{nf}} \right)^{\frac{2}{3}} k_B v$$

where, k_{bm} is the thermal conductivity value obtained through the modified Bridgman equation, ρ_{nf} is the density of nanofluid, and $M_{nf} = x_{bf} M_{bf} + x_p M_p$ is the molar mass of nanofluid. x_{bf} and x_p are the molar fractions of the base fluid and nanoparticle respectively whereas M_{bf} and M_p are the respective molar masses of the base fluid and nanoparticle.

WORKING PRINCIPLE



Ultrasound waves of known frequency are produced and its wavelength is measured. Then sound velocity in Liquid

$$v = \lambda \cdot f$$

After calculating velocity of sound in Liquid, one can calculate the thermal conductivity by the formula given by P W Bridgman.

DESCRIPTION

Thermal Conductivity Apparatus consists of following parts: Electronic Unit, Conductivity Cell-2MHz, Stability Cell 4MHz to increase settling time of the suspension, Temperature Controller Unit-To maintain temp. of liquids / nanofluids at desired temp from RT to 70°C.



Courtesy: V.Sesha Sai Kumar

SALIENT FEATURES

1. Works upto 70°C
2. Provided with temp variation Unit.
3. Errors less than 5%
4. User Friendly

Manufactured By:



MITTAL ENTERPRISES®

2151/T-7C, New Patel Nagar, New Delhi – 110008

Telefax: 011-25702784; Mobile: +91-9810681132, +91-9868532156

E-mail : info@mittalenterprises.com, mittalenterprises@bol.net.in,

Website : <http://www.mittalenterprises.com/>